

°CRYO Penguin

PERFORMANCE

FEATURES

- 1 Protocols for °CRYO Facials, Localised °CRYO, and Cellulite Reduction treatments
- 2 LED lights and precision lasers to indicate optimal distance, pinpoint treatment area, and alert user when critical skin temperature is reached
- 3 Ergonomic handle with remote control feature
- 4 Auto-balancer to maximize ease-of-use
- 5 10" Digital touchscreen with user-friendly interface
- 6 50 liter mobile nitrogen vessel that can be refilled at any time

CRYO SCIENCE sp. z o.o.
contact@cryosc.com

www.cryosc.com



LOCALIZED CRYOTHERAPY MACHINE

°CRYO Penguin™ is designed to deliver the safest and most effective localized cryotherapy treatments in the industry. With its ergonomic design, pre-programmed treatment protocols and built-in safety features, the device exceeds the highest standards in cryotherapy treatments.

TECHNICAL SPECIFICATIONS

Cooling Medium	Liquid Nitrogen (LN2)
Vessel Capacity	50 liters
Temperature of vapors of LN2 – smoothly regulated	Approx. -160°C
User Interface	Touch Screen 9.7”
Power	110/230V, 50/60 Hz
Intensity of Nitrogen vapour	Smoothly regulated

